

Glyphosate in Yellow Peas and Red Lentils Sample Preparation

1. Intended Use

For the detection of Glyphosate in ground yellow peas and red lentils.

2. Sensitivity

7.5 ppb in matrix

3. Materials and Reagents Required

Analytical balance

Weighing boats

Spatula

Microcentrifuge tubes, 1.5 mL or 2.0 mL

Microcentrifuge capable of 10,000 rpm

Micropipettes with disposable plastic tips

Glass vials – 4 mL, 20 mL with Teflon-lined caps

Blender/food processor or equivalent

Deionized water

Serological pipettes, 5 mL or 10 mL

Rotator and/or shaker

Vortex mixer

ABRAXIS® Glyphosate Sample Diluent (PN 500082)

ABRAXIS® Glyphosate Plate ELISA Kit (PN 500205)

4. Notes and Precautions

This procedure is intended for use with ground yellow peas and red lentils. Other matrices should be thoroughly validated before use with this procedure.

- It is highly recommended for the samples to be ground before extraction to produce accurate results. Ground the samples with a blender, food processor, or equivalent. Proceed to step 5.1.
- To minimize the potential for sample contamination due to carryover from the preparation of highly contaminated samples, the blender jar should be thoroughly washed and dried between samples.
- Due to the viscous nature of the sample extracts, if possible, the microtiter plate should be placed on a plate shaker or vortex mixer fitted with a micro-well plate holder adapter for the incubations with the antibody and conjugate solutions (Steps F.2 and F.3 in the ABRAXIS® Glyphosate ELISA user's guide). This will allow for the appropriate mixing of all reagents in the microtiter wells.
- Analysis should be performed with the ABRAXIS® Glyphosate Plate ELISA Kit as soon as possible after extraction. Samples should not sit more than one day in plastic microcentrifuge tubes before being diluted and analyzed.
- This procedure is for research use only. It is not intended for diagnostic procedures.

5. Procedure

5.1 Weigh 0.5 g of ground grain or flour samples to 20 mL glass vial.

5.2 Add 10 mL of deionized water to samples (1:20 dilution).

5.3 Vortex vigorously for 10 – 15 seconds and put samples on rotator or shaker at 40 rpm or 50% speed for 10 minutes.

5.4 Remove from rotator or shaker and allow the sample to settle for at least 2 minutes.

5.5 Transfer 1.5 to 2 mL of the supernatant to an appropriately labeled microcentrifuge vial.

5.6 Centrifuge for 5 minutes at ~8000 x g. Make sure the centrifuge is properly balanced.

- 5.7 Add 800 μ L of ABRAXIS[®] Glyphosate Sample Diluent to an appropriately labeled 4 mL glass vial. Add 200 μ L of the supernatant (from 5.6) to the ABRAXIS[®] Glyphosate Diluent in the vial (1:5 dilution). Vortex.
- 5.8 Derivatize the sample according to Section D Test Preparation in step 7 of *Derivatization of Standards, Control and Samples* instructions of the ABRAXIS[®] Glyphosate ELISA kit.
- 5.9 Perform assay as noted in Section F Assay Procedure instructions provided in the kit.

6. Evaluation of Results

The ELISA results must be multiplied by a factor of 100 to account for the necessary dilution. Samples showing a concentration lower than Standard 1 (0.075 ppb) should be reported as < 7.5 ppb of Glyphosate. Samples showing a higher concentration than Standard 5 (4.0 ppb) can be reported as > 400 ppb or diluted further and re-analyzed to obtain an accurate quantitative result.

7. For ordering or technical assistance contact

Gold Standard Diagnostics

124 Railroad Drive

Warminster, PA 18974

WEB: www.abraxiskits.com

Phone: (215) 357 3911

Fax: (215) 357 5232

Ordering: info.abraxaxis@us.goldstandarddiagnostics.com

Technical Support: support.abraxaxis@us.goldstandarddiagnostics.com

Date this Technical Bulletin is effective: 28FEB2022

Version: 01