



Contract Services and Manufacturing

Our scientific staff has decades of experience developing and manufacturing immunoassay products. This expertise allows us to provide customized research, development and production services that meet your need for high quality standards at competitive pricing. Below are examples of our key services and capabilities.

Services

- High-volume **microtiter plate coating**
- Antibody and antigen **production** and **purification**
- **Immunogold** reagents for rapid test system applications
- Bulk and custom **magnetic particles**
- **Receptor binding assays**, including bulk coated plates
- **ELISA** assay development
- **OEM manufacture** of kits, solutions, etc.

Capabilities

- Diversified test format (magnetic particles, microtiter plates, coated tubes, test strips, etc.)
- GLP documentation of product development and validation processes
- Quality manufacturing to your specifications
- Data evaluation options, from screening to fully quantitative configurations
- Design and validation of methods for sample matrices including foods, soils, water, etc.

Service Highlight: With a capacity of 2,000 microtiter plates per run and features such as liquid level sensing, bubble sensing and bar code printing, our high-volume microtiter plate coating services are completed quickly and to your exact specifications. Our technical expertise includes coating density optimization, stability maximization and the coating of surfaces with antibodies, receptors and other molecules.